


Special fume cupboards

Special application fume cupboard

Intended use

- Protective device for the user, tested in acc. with DIN 12924-2
- Suitable for open, thermal processes of breaking down chemicals with aggressive media such as e. g. sulphuric acid, perchloric acid, hydrofluoric acid or aqua regia
- The construction of the fume cupboard and the materials of the inner lining of the internal workspace determine which aggressive media the device can be used for
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- To prevent the formation of dangerous potentially explosive atmospheres in the internal workspace
- Protection from splashes of hazardous substances in the internal workspace
- Protection from flying particles, bodies or parts escaping from the internal workspace
- Fume cupboards constructed in acc. with DIN 12924-2, are normally not permitted for use with radioactive substances or microorganisms

Design




- 1 Sash with handle
- 2 Worktop
- 3 FAZ or AC control panel
- 4 Removable fascia panel
- 5 Extract air spigot integrated in fume-scrubber for harmful gases
- 6 Baffle
- 7 H-frame with push-in under-bench unit with support and service panels

Special fume cupboards

Special application fume cupboard

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 |
|---------------------------------------|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 |
| Depth [mm] | 900 | | |
| Height [mm] | 2700 | | |
| Clear width, internal workspace [mm] | 1150 | 1450 | 1750 |
| Clear height, internal workspace [mm] | 1125 | | |
| Working height [mm] | 900 | | |

| Weight | 1200 | 1500 | 1800 |
|--|-----------------|-------------|-------------|
| Without installations and fume-scrubber [kg] | Approx. 250 | Approx. 300 | Approx. 350 |
| Fume-scrubber without filling [kg] | 90 (type C 54) | | |
| | 100 (type C 90) | | |

Special fume cupboards

Special application fume cupboard

| Design characteristics | |
|--|---------------------------------------|
| Supporting construction | H-frame with push-in underbench units |
| Fume-scrubber | Optional |
| Extract manifold with condensate drain | Optional |
| Extract manifold with sprinkler | Optional |
| Neutralisation unit underbench unit | Optional |

| Electrics | |
|--------------------|------------------------------------|
| Electrical supply | External sockets in service panels |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | With take-off valves for vacuum, gases and/or waters and drip cup integrated in the worktop as an option |

| Ventilation technology | 1200 | 1500 | 1800 |
|---|---|----------------|----------------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 600 | 750 | 900 |
| Pressure loss, extract manifold with condensate drain [Pa] | 45/120 | 50/120 | 85/150 |
| Pressure loss, extract manifold with FAZ/AC [Pa] | FAZ 250/AC 300 | FAZ 300/AC 350 | FAZ 440/AC 500 |
| Pressure loss, fume cupboard with fume-scrubber [Pa] | 410/460 | 460/510 | 850/900 |
| Fume-scrubber Friatec | C 54 | | C 90 |
| Function display | FAZ | | |
| Airflow damper, constant | Airflow-Controller AC | | |
| Connection height [mm] for FAZ and AC with extract air spigot Ø 250 mm with fume-scrubber | 3140 | | |
| Connection height [mm] for FAZ with extract manifold Ø 250 mm without fume-scrubber | 2260 | | |
| Connection height [mm] for AC with extract manifold Ø 250 mm without fume-scrubber | 2490 | | |
| Underbench exhaust | As an option, depending on requirements and regulations | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175) and the maximum tracer gas values recommended by German Standard (BG Chemie).

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume cupboards with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|-----------------------------------|---|
| Internal lining including worktop | Stoneware (when sulphuric acid, aqua regia, perchloric acid are used) Polypropylene (when hydrofluoric acid is used) |


Special fume cupboards

Radio-isotope fume cupboard

Intended use

- Protective device for the user, tested in acc. with DIN 25466
- Extraction during work with radioactive substances if increased requirements for radiation protection apply
- Protection from incorporation, contamination and external radiation exposure
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- To prevent the formation of dangerous potentially explosive atmospheres in the internal workspace
- Protection from splashes of hazardous substances in the internal workspace
- Protection from flying particles, bodies or parts escaping from the internal workspace
- Fume cupboards constructed in acc. with DIN 25466 are normally not permitted for use with microorganisms
- Not suitable for openly breaking down chemicals

Design




- 1 Differential pressure gauge
- 2 Sash with handle
- 3 Worktop
- 4 FAZ or AC control panel
- 5 Removable fascia panel
- 6 Extract air spigot integrated in filter housing
- 7 Baffle with scaffold points
- 8 H-frame with push-in under-bench unit with support and service panels

Special fume cupboards

Radio-isotope fume cupboard

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 |
|---|-----------------|------|
| Width [mm] | 1200 | 1500 |
| Depth [mm] | 900 | |
| Height [mm] | 2700 | |
| Clear width, internal workspace [mm] | 1150 | 1450 |
| Clear height, internal workspace [mm] | 1053 | |
| Working height [mm] | 900 | |
| Filter housing, width x depth x height [mm] | 820 x 775 x 674 | |

| Weight | 1200 | 1500 |
|--|-------------|-------------|
| Without installations and lead insert [kg] | Approx. 250 | Approx. 300 |
| Filter housing [kg] | 90 | |

Special fume cupboards

Radio-isotope fume cupboard

| Design characteristics | |
|--|--|
| Supporting construction | Self-supporting underbench units or H-frame with push-in underbench units |
| Sash | One-piece |
| Number of devices for scaffold points, ø 12 to 13 mm | 6 |
| Filter, fume cupboard roof | Standard equipment: Filter F7 / particle filter H13 |
| Filter, lateral cabinet (max. 3 filter housings) | Filter housing, top: Particulate filter Filter housing, centre: Active charcoal filter Filter housing, bottom: Filter and particle filter |
| Differential pressure gauges | Display of the degree of saturation of the filters (not for active charcoal filter) |
| Lead insert | Optional |
| Waste disposal system for radio-isotope residual material in the underbench unit | Canister for liquid radio-isotope residual material as an option Collapsible boxes for solid radio-isotope residual material as an option Level indicator and/or opening in the worktop as an option |

| Electrics | |
|--------------------|------------------------------------|
| Electrical supply | External sockets in service panels |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | With take-off valves for vacuum and gases as an option |

| Ventilation technology | 1200 | 1500 |
|--|---|--------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 480 | 600 |
| Pressure loss, filter [Pa] ²⁾ | 25/200 | 30/235 |
| Pressure loss, particle filter [Pa] ²⁾ | 50/300 | 60/350 |
| Pressure loss, active charcoal filter [Pa] ²⁾ | 25/25 | 30/30 |
| Pressure loss, particulate filter [Pa] ²⁾ | 30/250 | 35/290 |
| Function display | FAZ | |
| Airflow damper, constant | Airflow-Controller AC | |
| Airflow damper, variable | Airflow-Controller AC | |
| Connection height [mm] for FAZ and AC with extract manifold Ø 250 mm | 3050 | |
| Underbench exhaust | As an option, depending on requirements and regulations | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175) and the maximum tracer gas values recommended by German Standard (BG Chemie).

²⁾ Pressure loss values refer to the states clean/contaminated.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume cupboards with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

In the case of fume cupboards with filters, the pressure loss of the integrated filter stages must be added to the pressure loss of the fume cupboard.

| Material/surface | |
|-----------------------------------|----------------------------------|
| Internal lining including worktop | Polypropylene Stainless steel |

Special fume cupboards

Radio-isotope fume cupboard

| Filter (filter in the filter cabinet or on the fume cupboard roof) | |
|---|--|
| Dimensions [mm] | 610 x 610 x 46 (+ 8 mm seal) |
| Pressure loss [Pa] at 1900 m ³ /h | 110 |
| Design characteristics | Filter element (fine particle filter); filter class EN 779: F7 Frame made of multilayer board with grip and type label on the 610-mm side; PU seal on the dust-laden air side |
| Use | Fine particle filter for particle adsorption, e.g.: Oil smoke and agglomerated soot, tobacco smoke, metal oxide smoke Average efficiency (Em) 80–90% |

| Particle filter (filter in the filter cabinet or on the fume cupboard roof) | |
|--|---|
| Dimensions [mm] | 610 x 610 x 292 (+ 7 mm seal) |
| Pressure loss [Pa] at 2435 m ³ /h | 250 |
| Design characteristics | Particle filter element type: Hepa H13; efficiency: MPPS Frame made of multilayer board with grip and type label on the 610-mm side; PU tight seat seal on the clean air side; filter medium flush on the clean air side |
| Use | Particle filter for the adsorption of particles up to H13; particle adsorption 99.95 %; transmittance 0.05% |

| Active charcoal filter (filter in the filter cabinet) | |
|--|---|
| Dimensions [mm] | 610 x 610 x 292 (+ 7 mm seal) |
| Pressure loss [Pa] at 600 m ³ /h | 9 |
| Design characteristics | Activated charcoal cell 7C for 16 x activated charcoal cartridges Frame galvanised sheet metal; 2 x U handle and type plate on the 610-mm side; PU tight seat seal on the clean air side |
| Use | Standard impregnation: for all common radioactive materials, radioactive iodine compounds, radioactive iodomethane, radioactive gases. (A filter with filters class F7 in acc. with EN 779 is recommended.) |

| Particulate filter (filter in the filter cabinet) | |
|--|--|
| Dimensions [mm] | 610 x 610 x 292 (+ 7 mm seal) |
| Pressure loss [Pa] at 1965 m ³ /h | 125 |
| Design characteristics | Particulate or Micretain filter element type: Hepa H11 in acc. with EN 1822 Frame made of multilayer board with grip and type label on the 610-mm side; PU tight seat seal on the clean air side; filter medium flush on the clean air side |
| Use | Particle filter for the adsorption of particles up to H11; particle adsorption 95 %; transmittance 5%; to be installed after active charcoal filters to bind the charcoal dust contamination from the charcoal filter. |


Special fume cupboards

Filter fume cupboard

Intended use

- Before the extract air from the internal workspace is released into the environment, it is cleaned by a filter unit

Design




- 1 Differential pressure gauge
- 2 Sash with handle and horizontal sashes
- 3 Worktop
- 4 FAZ or AC control panel
- 5 Removable fascia panel
- 6 Extract air spigot
- 7 Baffle with scaffold points
- 8 H-frame with push-in under-bench unit with support and service panels

Special fume cupboards

Filter fume cupboard

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 |
|---|-----------------|------|------|
| Width [mm] | 1200 | 1500 | 1800 |
| Depth [mm] | 900 | | |
| Height [mm] | 2700 | | |
| Clear width, internal workspace [mm] | 1150 | 1450 | 1750 |
| Clear height, internal workspace [mm] | 1060 | | |
| Working height [mm] | 900 | | |
| Filter housing, width x depth x height [mm] | 820 x 775 x 674 | | |

| Weight | 1200 | 1500 | 1800 |
|---|-------------|-------------|-------------|
| Filter fume cupboard without installations [kg] | Approx. 270 | Approx. 320 | Approx. 370 |
| Filter housing [kg] | 90 | | |

Special fume cupboards

Filter fume cupboard

| Design characteristics | 1200 | 1500 | 1800 |
|---|--|------|---------------------|
| Supporting construction | H-frame with push-in underbench units | | |
| Sash | 2 horizontal sashes | | 3 horizontal sashes |
| Glass pane in the side wall | Possible on the left and/or right side of the fume cupboard; not with stoneware internal lining | | |
| Number of devices for scaffold points, ø 12 to 13 mm | 6 | | 8 |
| Material lock | Possible on the left and/or right side of the fume cupboard | | |
| Filter, fume cupboard roof | Standard equipment: Filter F7 / particle filter H13 | | |
| Differential pressure gauges | Display of the degree of saturation of the filters | | |

| Electrics | |
|--------------------|------------------------------------|
| Electrical supply | External sockets in service panels |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | With take-off valves for vacuum, gases and/or waters and drip cup integrated in the worktop as an option |

| Ventilation technology | 1200 | 1500 | 1800 |
|---|---|--------|---------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 480 | 600 | 720 |
| Pressure loss, filter [Pa] ²⁾ | 35/200 | 45/235 | 65/290 |
| Pressure loss, particle filter [Pa] ²⁾ | 70/300 | 95/365 | 130/430 |
| Pressure loss, active charcoal filter [Pa] ²⁾ | 35/25 | 45/30 | 65/35 |
| Function display | FAZ | | |
| Airflow damper, constant | Airflow-Controller AC | | |
| Airflow damper, variable | Airflow-Controller AC | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | |
| Connection height [mm] for FAZ and AC with extract air spigot ø 250 mm | 3050 | | |
| Underbench exhaust | As an option, depending on requirements and regulations | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175) and the maximum tracer gas values recommended by German Standard (BG Chemie).

²⁾ Pressure loss values refer to the states clean/contaminated.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume cupboards with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

In the case of fume cupboards with filters, the pressure loss of the integrated filter stages must be added to the pressure loss of the fume cupboard.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Epoxy Stainless steel |
| Internal lining | Melamine resin facing Solid grade laminate |

Special fume cupboards

Filter fume cupboard

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Epoxy Stainless steel |
| Internal lining | Melamine resin facing Solid (grade) laminate |

| Filter | |
|--|--|
| Dimensions [mm] | 610 x 610 x 46 (+ 8 mm seal) |
| Pressure loss [Pa] at 1900 m ³ /h | 110 |
| Design characteristics | Filter element (fine particle filter); filter class EN 779: F7 Frame made of multilayer board with grip and type label on the 610-mm side; PU seal on the dust-laden air side |
| Use | Fine particle filter for particle adsorption, e.g.: Oil smoke and agglomerated soot, tobacco smoke, metal oxide smoke Average efficiency (Em) 80–90% |

| Particle filter | |
|--|---|
| Dimensions [mm] | 610 x 610 x 292 (+ 7 mm seal) |
| Pressure loss [Pa] at 2435 m ³ /h | 250 |
| Design characteristics | Particle filter element type: Hepa H13; efficiency: MPPS Frame made of multilayer board with grip and type label on the 610-mm side; PU tight seat seal on the clean air side; filter medium flush on the clean air side |
| Use | Particle filter for the adsorption of particles up to H13; particle adsorption 99.95 %; transmittance 0.05% |