


## Top-mounted cabinets

### Top-mounted cabinet

#### Intended use

- For storing equipment and chemicals in acc. with EN 14727
- Only suitable as a permanently installed top part on the following Waldner cabinets: Laboratory cabinet, pull-out cabinet, emergency cabinet and acids and alkalis cabinet
- Not permitted for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not permitted for storing acids and alkalis

#### Design



#### Technical data

Dimensions				
Width [mm]	450	600	900	1200
Depth [mm]			350	550
Height [mm]			610	760

Load bearing capacity	
Per shelf [kg]	30

Design characteristics	
Construction	With integrated rail for securing a ladder For laboratory cabinets with or without extract air spigot 1 shelf, height-adjustable Hinged doors
Handle	U handle SCALA U handle, stainless steel
Hook ladder	Optional
Closing	Optional