


Pull-out cabinets

Pull-out cabinet

Intended use

- For storing liquid or solid substances in suitable containers in acc. with EN 14727
- Not permitted for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not permitted for storing acids and alkalis

Design



- 1 Wire basket with tray
- 2 Pull-out

Pull-out cabinets

Pull-out cabinet

Technical data

Dimensions	
Width [mm]	600 900
Depth [mm]	550
Overall height [mm]	2090
Height, plinth [mm]	110
Tray, width x depth x height [mm]	240 x 425 x 40

Load bearing capacity	
Per drawer [kg]	120
Per tray [kg]	10

Design characteristics	
Construction	5 wire baskets with trays for each drawer, height-adjustable Fastened to the wall 4 height-adjustable feet Drawer doors with drawers accessible from both sides
Handle	U handle SCALA U handle, stainless steel
Soft stop for drawers	Optional
Compartment partitioning	Optional
Extract air spigot	Optional
Closing	Optional

Material	
Tray	Polyethylene